

For Research Use Only

Galectin-3 Monoclonal antibody

Catalog Number: 60207-2-Ig **1 Publications**



Basic Information

Catalog Number: 60207-2-Ig	GenBank Accession Number: BC001120	Purification Method: Protein G purification
Size: N/A	GeneID (NCBI): 3958	CloneNo.: 3H10B6
Source: Mouse	UNIPROT ID: P17931	
Isotype: IgG1	Full Name: lectin, galactoside-binding, soluble, 3	
Immunogen Catalog Number: AG6753	Calculated MW: 26 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human, mouse

Cited Species:
mouse

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data