

For Research Use Only

phospho(409/410)-TDP43 Monoclonal antibody



Catalog Number: 60218-1-Ig

Basic Information

Catalog Number:

60218-1-Ig

Size:

600 µg/ml

Source:

Mouse

Isotype:

IgM

GenBank Accession Number:

NM_007375

GeneID (NCBI):

23435

UNIPROT ID:

Q13148

Full Name:

TAR DNA binding protein

Calculated MW:

43 kDa

Purification Method:

Caprylic acid/ammonium sulfate precipitation

CloneNo.:

1E9A12

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data