

For Research Use Only

CD44 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60224-2-PBS



Basic Information

Catalog Number: 60224-2-PBS	GenBank Accession Number: BC004372	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 960	CloneNo.: 1G7G4
Source: Mouse	ENSEMBL Gene ID: ENSG00000026508	
Isotype: IgG1	UNIPROT ID: P16070	
Immunogen Catalog Number: AG7633	Full Name: CD44 molecule (Indian blood group)	
	Calculated MW: 742 aa, 82 kDa	

Applications

Tested Applications:
IF/ICC, Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50724-1, CD44 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60224-2-PBS. Detection antibody: 60224-3-PBS. Standard:Ag7633. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50724-3, CD44 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60224-2-PBS. Detection antibody: 60224-4-PBS. Standard:Ag7633. Range: 0.098-100 ng/mL