

For Research Use Only

PFKFB3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60241-2-PBS



Basic Information

Catalog Number:

60241-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4737

GenBank Accession Number:

BC040482

GeneID (NCBI):

5209

UNIPROT ID:

Q16875

Full Name:

6-phosphofructo-2-kinase/fructose-
2,6-biphosphatase 3

Calculated MW:

520 aa, 60 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3E12A10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51513-1, PFKFB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60241-2-PBS. Detection antibody: 60241-1-PBS.



Cytometric bead array standard curve of MP51513-1, PFKFB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60241-2-PBS. Detection antibody: 60241-1-PBS. Standard:Ag4737. Range: 3.125-400 ng/mL.