

For Research Use Only

OCT4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60242-6-PBS



Basic Information

Catalog Number:

60242-6-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1794

GenBank Accession Number:

BC020712

GeneID (NCBI):

5460

UNIPROT ID:

Q01860

Full Name:

POU class 5 homeobox 1

Calculated MW:

39 kDa

Purification Method:

Protein G purification

CloneNo.:

5H9G12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50231-3, OCT4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60242-5-PBS. Detection antibody: 60242-6-PBS. Standard:Ag1794. Range: 0.391-25 ng/mL