

For Research Use Only

CDX2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60243-3-PBS



Basic Information

Catalog Number:

60243-3-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG17644

GenBank Accession Number:

BC014461

GeneID (NCBI):

1045

UNIPROT ID:

Q99626

Full Name:

caudal type homeobox 2

Calculated MW:

313 aa, 34 kDa

Purification Method:

Protein A purification

CloneNo.:

3A1G8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50121-1, CDX2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60243-1-PBS. Detection antibody: 60243-3-PBS. Standard:Ag17644. Range: 0.098-100 ng/mL.