

For Research Use Only

# KRT15 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 60247-2-PBS

## Basic Information

Catalog Number: 60247-2-PBS	GenBank Accession Number: BC002641	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 3866	CloneNo.: 4B7B6
Source: Mouse	UNIPROT ID: P19012	
Isotype: IgG2b	Full Name: keratin 15	
Immunogen Catalog Number: AG0185	Calculated MW: 49 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Sandwich ELISA standard curve ofMP50002-1, KRT15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60247-1-PBS. Detection antibody: HRP-conjugated 60247-2-PBS. Standard: Ag0185. Range: 0.313-20 ng/mL



Cytometric bead array standard curve ofMP50002-1, KRT15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60247-1-PBS. Detection antibody: 60247-2-PBS. Standard: Ag0185. Range: 1.563-100 ng/mL