

For Research Use Only

SND1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60265-2-PBS



Basic Information

Catalog Number: 60265-2-PBS	GenBank Accession Number: BC017180	Purification Method: Protein A purification
Source: Mouse	GeneID (NCBI): 27044	CloneNo.: 2D5F6
Isotype: IgG2b	UNIPROT ID: Q7KZF4	
Immunogen Catalog Number: AG1200	Full Name: staphylococcal nuclease and tudor domain containing 1	
	Calculated MW: 101 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51825-1, SND1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60265-1-PBS. Detection antibody: 60265-2-PBS.



Cytometric bead array standard curve of MP51825-1, SND1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60265-1-PBS. Detection antibody: 60265-2-PBS. Standard:Ag1200. Range: 0.156-20 ng/mL.