

For Research Use Only

Podoplanin Monoclonal antibody

Catalog Number: 60275-1-Ig



Basic Information

Catalog Number: 60275-1-Ig	GenBank Accession Number: BC022812	Purification Method: Protein G purification
Size: 2000 µg/ml	GeneID (NCBI): 10630	CloneNo.: 1F2A9
Source: Mouse	UNIPROT ID: Q86YL7	
Isotype: IgG1	Full Name: podoplanin	
Immunogen Catalog Number: AG2201	Calculated MW: 17 kDa, 25 kDa	

Applications

Tested Applications:
IHC, WB, ELISA

Species Specificity:
human, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data