

For Research Use Only

Alpha Synuclein Monoclonal antibody, PBS Only (Capture)



Catalog Number: 60285-1-PBS

Basic Information

Catalog Number:

60285-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1285

GenBank Accession Number:

BC013293

GeneID (NCBI):

6622

UNIPROT ID:

P37840

Full Name:

synuclein, alpha (non A4 component of amyloid precursor)

Calculated MW:

14 kDa

Observed MW:

18-20 kDa

Purification Method:

Protein A purification

CloneNo.:

1B8A5

Applications

Tested Applications:

ELISA

Species Specificity:

human, rat

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS only

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP50043-1, α -Synuclein Monoclonal Matched Antibody Pair - BSA and Azide free, measured by sandwich ELISA. Capture antibody: 60285-1-PBS. Detection antibody: 66412-1-PBS (HRP conjugated). Standard: Ag1285. Range: 0.313-20 ng/mL