

For Research Use Only

Alpha Synuclein Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60285-1-PBS



Basic Information

Catalog Number:

60285-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1285

GenBank Accession Number:

BC013293

GeneID (NCBI):

6622

UNIPROT ID:

P37840

Full Name:

synuclein, alpha (non A4 component
of amyloid precursor)

Calculated MW:

14 kDa

Observed MW:

18-20 kDa

Purification Method:

Protein G purification

CloneNo.:

1B8A5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50043-1, α-Synuclein Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60285-1-PBS. Detection antibody: 66412-1-PBS. Standard: Ag1285. Range: 0.313-20 ng/mL.