

For Research Use Only

CD45 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60287-2-PBS



Basic Information

Catalog Number:

60287-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG13712

GenBank Accession Number:

BC014239

GeneID (NCBI):

5788

ENSEMBL Gene ID:

ENSG00000081237

UNIPROT ID:

P08575

Full Name:

protein tyrosine phosphatase,
receptor type, C

Calculated MW:

1304 aa, 147 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

4F3B7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50221-1, CD45 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60287-2-PBS. Detection antibody: 60287-3-PBS. Standard:Ag13712. Range: 1.563-100 ng/mL.