

For Research Use Only

U2AF35 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60289-2-PBS



Basic Information

Catalog Number:

60289-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0399

GenBank Accession Number:

BC001923

GeneID (NCBI):

7307

UNIPROT ID:

Q01081

Full Name:

U2 small nuclear RNA auxiliary factor

1

Calculated MW:

28 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

6D3D12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51496-1, U2AF35 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60289-2-PBS. Detection antibody: 60289-1-PBS.



Cytometric bead array standard curve of MP51496-1, U2AF35 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60289-2-PBS. Detection antibody: 60289-1-PBS. Standard:Ag0399. Range: 3.125-400 ng/mL