

For Research Use Only

BSA Monoclonal antibody

Catalog Number: 60302-1-Ig **1 Publications**



Basic Information

Catalog Number:

60302-1-Ig

Size:

1800 µg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

GeneID (NCBI):

280717

Full Name:

bovine serum albumin

Observed MW:

66 kDa

Purification Method:

Protein G purification

CloneNo.:

7C6B5

Applications

Tested Applications:

ELISA

Species Specificity:

bovine

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

