

For Research Use Only

Caspase 1/p20/p10 Monoclonal antibody

Catalog Number: 60313-1-Ig



Basic Information

Catalog Number:

60313-1-Ig

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11143

GenBank Accession Number:

BC062327

GeneID (NCBI):

834

UNIPROT ID:

P29466

Full Name:

caspase 1, apoptosis-related cysteine peptidase (interleukin 1, beta, convertase)

Calculated MW:

404 aa, 45 kDa

Purification Method:

Protein G purification

CloneNo.:

5H2B4

Applications

Tested Applications:

ELISA

Species Specificity:

human, rat

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

