

For Research Use Only

hCG Beta Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60334-2-PBS



Basic Information

Catalog Number:

60334-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG2191

GenBank Accession Number:

BC022796

GeneID (NCBI):

1082

UNIPROT ID:

PODN86

Full Name:

CG beta

Calculated MW:

165 aa, 18 kDa

Purification Method:

Protein A purification

CloneNo.:

2F3D4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50292-1, CGB,hCG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60334-2-PBS. Detection antibody: 60334-1-PBS. Standard:Ag2191. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50292-2, CGB,hCG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60334-2-PBS. Detection antibody: 60334-3-PBS. Standard:Ag2191. Range: 0.098-100 ng/mL.