

For Research Use Only

BSAP,PAX5 Monoclonal antibody, PBS Only (Detector)



Catalog Number:60349-3-PBS

Basic Information

Catalog Number:

60349-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG16023

GenBank Accession Number:

BC156927

GeneID (NCBI):

5079

UNIPROT ID:

Q02548

Full Name:

paired box 5

Calculated MW:

391 aa, 42 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

5E1E7

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50228-1, BSAP,PAX5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60349-2-PBS. Detection antibody: 60349-3-PBS. Standard:Ag16023. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP50228-2, BSAP,PAX5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60349-1-PBS. Detection antibody: 60349-3-PBS. Standard:Ag16023. Range: 3.125-100 ng/mL