

For Research Use Only

LMOD1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60414-2-PBS



Basic Information

Catalog Number:

60414-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG7568

GenBank Accession Number:

BC001755

GeneID (NCBI):

25802

UNIPROT ID:

P29536

Full Name:

leiomodin 1 (smooth muscle)

Calculated MW:

67 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1C7D3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51375-1, LMOD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60414-2-PBS. Detection antibody: 60414-1-PBS.



Cytometric bead array standard curve of MP51375-1, LMOD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60414-2-PBS. Detection antibody: 60414-1-PBS. Standard: Ag7568. Range: 0.195-200 ng/mL