

For Research Use Only

MTIF3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60416-1-PBS



Basic Information

Catalog Number: 60416-1-PBS	GenBank Accession Number: BC046166	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 219402	CloneNo.: 2F8B4
Isotype: IgG1	UNIPROT ID: Q9H2K0	
Immunogen Catalog Number: AG5643	Full Name: mitochondrial translational initiation factor 3	
	Calculated MW: 32 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50533-1, MTIF3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60416-1-PBS. Detection antibody: 60416-2-PBS. Standard:Ag5643. Range: 0.781-100 ng/mL.



Biolayer interferometry (BLI) kinetic assays of 60416-1-Ig against Human MTIF3 were performed. The affinity constant is below 1 pM.