

For Research Use Only

SLC10A3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60419-1-PBS



Basic Information

Catalog Number:

60419-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG14038

GenBank Accession Number:

BC004966

GeneID (NCBI):

8273

UNIPROT ID:

P09131

Full Name:

solute carrier family 10 (sodium/bile acid cotransporter family), member 3

Calculated MW:

477 aa, 50 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H1A11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50537-1, SLC10A3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60419-1-PBS. Detection antibody: 60419-2-PBS. Standard:Ag14038. Range: 0.781-12.5 ng/mL