

For Research Use Only

MYH1 Monoclonal antibody, PBS Only

Catalog Number: 60421-2-PBS



Basic Information

Catalog Number:

60421-2-PBS

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG27992

GenBank Accession Number:

BC114545

GeneID (NCBI):

4619

UNIPROT ID:

P12882

Full Name:

myosin, heavy chain 1, skeletal muscle, adult

Calculated MW:

1939 aa, 223 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3A5A4

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

