

For Research Use Only

CLTA Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60427-4-PBS



Basic Information

Catalog Number:

60427-4-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1299

GenBank Accession Number:

BC019287

GeneID (NCBI):

1211

UNIPROT ID:

P09496

Full Name:

clathrin, light chain (Lca)

Calculated MW:

27 kDa

Purification Method:

Protein G purification

CloneNo.:

2F9D4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50560-4, CLTA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60427-3-PBS. Detection antibody: 60427-4-PBS. Standard:Ag1299. Range: 1.563-100 ng/mL.



Biolayer interferometry (BLI) kinetic assays of 60427-4-Ig against Human CLTA were performed. The affinity constant is 3.29 nM.