

For Research Use Only

# CPVL Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60429-2-PBS



## Basic Information

Catalog Number:

60429-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG36547

GenBank Accession Number:

BC016838

GeneID (NCBI):

54504

UNIPROT ID:

Q9H3G5

Full Name:

carboxypeptidase, vitellogenic-like

Calculated MW:

476 aa, 54 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E6E8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


# Selected Validation Data



Cytometric bead array standard curve of MP50562-1, CPVL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60429-1-PBS. Detection antibody: 60429-2-PBS. Standard:Ag36547. Range: 3.125-100 ng/mL



Cytometric bead array standard curve of MP50562-2, CPVL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60429-3-PBS. Detection antibody: 60429-2-PBS. Standard:Ag36547. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP50562-4, CPVL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60429-2-PBS. Detection antibody: 60429-6-PBS. Standard:Ag36547. Range: 0.781-100 ng/mL