

For Research Use Only

CPVL Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60429-6-PBS



Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
60429-6-PBS	BC016838	Protein A Magarose purification
Size:	GeneID (NCBI):	CloneNo.:
1 mg/ml	54504	2D9H11
Source:	UNIPROT ID:	
Mouse	Q9H3G5	
Isotype:	Full Name:	
IgG2a	carboxypeptidase, vitellogenic-like	
Immunogen Catalog Number:	Calculated MW:	
AG36547	476 aa, 54 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50562-4, CPVL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60429-2-PBS. Detection antibody: 60429-6-PBS. Standard:Ag36547. Range: 0.781-100 ng/mL.