

For Research Use Only

# CYP11B1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60435-2-PBS



## Basic Information

Catalog Number:

60435-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG17158

GenBank Accession Number:

BC096285

GeneID (NCBI):

1584

UNIPROT ID:

P15538

Full Name:

cytochrome P450, family 11, subfamily B, polypeptide 1

Calculated MW:

574 aa, 64 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1C1B10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50582-1, CYP11B1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60435-1-PBS. Detection antibody: 60435-2-PBS. Standard: Ag17158. Range: 0.098-50 ng/mL