

For Research Use Only

CYP11B1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60435-5-PBS



Basic Information

Catalog Number:

60435-5-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17158

GenBank Accession Number:

BC096285

GeneID (NCBI):

1584

UNIPROT ID:

P15538

Full Name:

cytochrome P450, family 11, subfamily B, polypeptide 1

Calculated MW:

574 aa, 64 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G2C2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50582-3, CYP11B1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60435-3-PBS. Detection antibody: 60435-5-PBS. Standard:Ag17158. Range: 0.098-50 ng/mL



Cytometric bead array standard curve of MP50582-4, CYP11B1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60435-1-PBS. Detection antibody: 60435-5-PBS. Standard:Ag17158. Range: 0.098-50 ng/mL