

For Research Use Only

NSUN4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60448-2-PBS



Basic Information

Catalog Number:

60448-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG9417

GenBank Accession Number:

BC014441

GeneID (NCBI):

387338

UNIPROT ID:

Q96CB9

Full Name:

NOL1/NOP2/Sun domain family, member 4

Calculated MW:

384 aa, 43 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1E9B5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50604-1, NSUN4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60448-1-PBS. Detection antibody: 60448-2-PBS. Standard:Ag9417. Range: 0.195-100 ng/mL.