

For Research Use Only

PLEKHJ1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60453-5-PBS



Basic Information

Catalog Number:

60453-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG13729

GenBank Accession Number:

BC003084

GeneID (NCBI):

55111

UNIPROT ID:

Q9NW61

Full Name:

pleckstrin homology domain
containing, family J member 1

Calculated MW:

149 aa, 18 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2D4E5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50615-3, PLEKHJ1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60453-5-PBS. Detection antibody: 60453-2-PBS. Standard:Ag13729. Range: 3.125-50 ng/mL