

For Research Use Only

# B3GAT1 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60455-2-PBS



## Basic Information

Catalog Number:

60455-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

NM\_018644.3

GeneID (NCBI):

27087

UNIPROT ID:

Q9P2W7-1

Full Name:

beta-1,3-glucuronyltransferase 1  
(glucuronosyltransferase P)

Calculated MW:

38 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1C8D12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50618-1, B3GAT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60455-1-PBS. Detection antibody: 60455-2-PBS. Standard:Eg0793. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50618-3, B3GAT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60455-2-PBS. Detection antibody: 60455-3-PBS. Standard:Eg0793. Range: 0.195-100 ng/mL