

For Research Use Only

CPN1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60459-1-PBS



Basic Information

Catalog Number:

60459-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG36462

GenBank Accession Number:

BC027897

GeneID (NCBI):

1369

UNIPROT ID:

P15169

Full Name:

carboxypeptidase N, polypeptide 1

Calculated MW:

458 aa, 52 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3C3B2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50626-1, CPN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60459-1-PBS. Detection antibody: 60459-2-PBS. Standard:Ag36462. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP50626-2, CPN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60459-1-PBS. Detection antibody: 60459-3-PBS. Standard:Ag36462. Range: 1.563-100 ng/mL