

For Research Use Only

MCEE Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60461-1-PBS



Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
60461-1-PBS	BC020825	Protein G Magarose purification
Size:	GeneID (NCBI):	CloneNo.:
1 mg/ml	84693	1G2C2
Source:	UNIPROT ID:	
Mouse	Q96PE7	
Isotype:	Full Name:	
IgG1	methylmalonyl CoA epimerase	
Immunogen Catalog Number:	Calculated MW:	
AG9665	176 aa, 19 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50633-1, MCEE Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60461-1-PBS. Detection antibody: 60461-2-PBS. Standard:Ag9665. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50633-2, MCEE Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60461-1-PBS. Detection antibody: 60461-3-PBS. Standard:Ag9665. Range: 0.391-100 ng/mL