

For Research Use Only

RTN4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60462-1-PBS



Basic Information

Catalog Number:

60462-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33806

GenBank Accession Number:

BC007109

GeneID (NCBI):

57142

UNIPROT ID:

Q9NQC3

Full Name:

reticulon 4

Calculated MW:

130 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2E10G1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50635-1, RTN4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60462-1-PBS. Detection antibody: 60462-2-PBS. Standard:Ag33806. Range: 0.098-6.25 ng/mL



Cytometric bead array standard curve of MP50635-2, RTN4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60462-1-PBS. Detection antibody: 60462-3-PBS. Standard:Ag33806. Range: 0.098-100 ng/mL