

For Research Use Only

CCNF Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60470-1-PBS



Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
60470-1-PBS	BC012349	Protein G Magarose purification
Size:	GeneID (NCBI):	CloneNo.:
1 mg/ml	899	1G10C6
Source:	UNIPROT ID:	
Mouse	P41002	
Isotype:	Full Name:	
IgG1	cyclin F	
Immunogen Catalog Number:	Calculated MW:	
AG9322	786 aa, 88 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50640-1, CCNF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60470-1-PBS. Detection antibody: 60470-2-PBS. Standard:Ag9322. Range: 0.195-12.5 ng/mL