

For Research Use Only

TMIGD1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60472-4-PBS



Basic Information

Catalog Number:

60472-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG19612

GenBank Accession Number:

BC137201

GeneID (NCBI):

388364

UNIPROT ID:

Q6UXZ0

Full Name:

transmembrane and immunoglobulin
domain containing 1

Calculated MW:

262 aa, 29 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2F8A7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50643-2, TMIGD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60472-3-PBS. Detection antibody: 60472-4-PBS. Standard:Ag19612. Range: 1.563-100 ng/mL.