

For Research Use Only

MRPL48 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60474-1-PBS



Basic Information

Catalog Number:

60474-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6819

GenBank Accession Number:

BC009764

GeneID (NCBI):

51642

UNIPROT ID:

Q96GC5

Full Name:

mitochondrial ribosomal protein L48

Calculated MW:

24 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A4C5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50648-1, MRPL48 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60474-1-PBS. Detection antibody: 60474-2-PBS. Standard:Ag6819. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50648-3, MRPL48 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60474-1-PBS. Detection antibody: 60474-5-PBS. Standard:Ag6819. Range: 0.781-100 ng/mL