

For Research Use Only

# PRAS40 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60477-1-PBS



## Basic Information

Catalog Number:

60477-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC051844

GeneID (NCBI):

84335

UNIPROT ID:

Q96B36

Full Name:

AKT1 substrate 1 (proline-rich)

Calculated MW:

27 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F2A9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50656-1, PRAS40 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60477-1-PBS. Detection antibody: 60478-1-PBS. Standard:Ag15355. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50656-2, PRAS40 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60477-1-PBS. Detection antibody: 60478-2-PBS. Standard:Ag15355. Range: 1.563-100 ng/mL.