

For Research Use Only

# MARS2 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60484-1-PBS



## Basic Information

Catalog Number:

60484-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG18596

GenBank Accession Number:

BC126294

GeneID (NCBI):

92935

UNIPROT ID:

Q96GW9

Full Name:

methionyl-tRNA synthetase 2,  
mitochondrial

Calculated MW:

593 aa, 67 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1A10F1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50670-1, MARS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60484-1-PBS. Detection antibody: 60484-2-PBS. Standard:Ag18596. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50670-2, MARS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60484-3-PBS. Detection antibody: 60484-1-PBS. Standard:Ag18596. Range: 1.563-100 ng/mL.