

For Research Use Only

# XPR1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60489-4-PBS



## Basic Information

Catalog Number:

60489-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5731

GenBank Accession Number:

BC041142

GeneID (NCBI):

9213

UNIPROT ID:

Q9UBH6

Full Name:

xenotropic and polytropic retrovirus receptor

Calculated MW:

696 aa, 82 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3D3C9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50679-2, XPR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60489-3-PBS. Detection antibody: 60489-4-PBS. Standard:Ag5731. Range: 1.563-200 ng/mL.