

For Research Use Only

RDH13 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60491-3-PBS



Basic Information

Catalog Number:

60491-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9107

GenBank Accession Number:

BC009881

GeneID (NCBI):

112724

UNIPROT ID:

Q8NBN7

Full Name:

retinol dehydrogenase 13 (all-trans/9-cis)

Calculated MW:

331 aa, 36 kDa

Purification Method:

Protein G purification

CloneNo.:

2A9C6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50684-2, RDH13 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60491-3-PBS. Detection antibody: 60491-2-PBS. Standard:Ag9107. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50684-3, RDH13 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60491-1-PBS. Detection antibody: 60491-3-PBS. Standard:Ag9107. Range: 0.098-100 ng/mL