

For Research Use Only

# MTHFD2L Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60504-2-PBS



## Basic Information

Catalog Number:

60504-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11470

GenBank Accession Number:

BC032771

GeneID (NCBI):

441024

UNIPROT ID:

Q9H903

Full Name:

methylenetetrahydrofolate  
dehydrogenase (NADP+ dependent) 2-  
like

Calculated MW:

146 aa, 16 kDa

Purification Method:

Protein G purification

CloneNo.:

1C1D7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50702-1, MTHFD2L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60504-1-PBS. Detection antibody: 60504-2-PBS. Standard:Ag11470. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50702-3, MTHFD2L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60504-3-PBS. Detection antibody: 60504-2-PBS. Standard:Ag11470. Range: 1.563-100 ng/mL