

For Research Use Only

MTHFD2L Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60504-3-PBS



Basic Information

Catalog Number:

60504-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11470

GenBank Accession Number:

BC032771

GeneID (NCBI):

441024

UNIPROT ID:

Q9H903

Full Name:

methylenetetrahydrofolate
dehydrogenase (NADP+ dependent) 2-
like

Calculated MW:

146 aa, 16 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A6B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50702-2, MTHFD2L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60504-1-PBS. Detection antibody: 60504-3-PBS. Standard:Ag11470. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP50702-3, MTHFD2L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60504-3-PBS. Detection antibody: 60504-2-PBS. Standard:Ag11470. Range: 1.563-100 ng/mL