

For Research Use Only

# PLD6 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60508-1-PBS



## Basic Information

Catalog Number:	60508-1-PBS	GenBank Accession Number:	BC031263	Purification Method:	Protein G purification
Size:	1 mg/ml	GeneID (NCBI):	201164	CloneNo.:	2E10G7
Source:	Mouse	UNIPROT ID:	Q8N2A8		
Isotype:	IgG1	Full Name:	phospholipase D family, member 6		
Immunogen Catalog Number:	AG33159	Calculated MW:	252 aa, 28 kDa		

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP50709-1, PLD6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60508-1-PBS. Detection antibody: 60508-2-PBS. Standard:Ag33159. Range: 0.098-100 ng/mL.