

For Research Use Only

TMEM63A Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60509-2-PBS



Basic Information

Catalog Number:

60509-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG16911

GenBank Accession Number:

BC030245

GeneID (NCBI):

9725

UNIPROT ID:

O94886

Full Name:

transmembrane protein 63A

Calculated MW:

807 aa, 92 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2D11H10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50711-1, TMEM63A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60509-1-PBS. Detection antibody: 60509-2-PBS. Standard:Ag16911. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP50711-2, TMEM63A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60509-3-PBS. Detection antibody: 60509-2-PBS. Standard:Ag16911. Range: 0.391-100 ng/mL