

For Research Use Only

# TCL1A Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60511-1-PBS



## Basic Information

Catalog Number: 60511-1-PBS	GenBank Accession Number: BC005831	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 8115	CloneNo.: 4G10G10
Source: Mouse	UNIPROT ID: P56279	
Isotype: IgG1	Full Name: T-cell leukemia/lymphoma 1A	
Immunogen Catalog Number: AG0786	Calculated MW: 13 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50713-1, TCL1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60511-1-PBS. Detection antibody: 60511-2-PBS. Standard:Ag0786. Range: 3.125-100 ng/mL.