

For Research Use Only

NTF3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60520-1-PBS



Basic Information

Catalog Number:

60520-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG35836

GenBank Accession Number:

BC069773

GeneID (NCBI):

4908

UNIPROT ID:

P20783

Full Name:

neurotrophin 3

Calculated MW:

257 aa, 29 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3G12F11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50760-1, NTF3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66933-2-PBS. Detection antibody: 60520-1-PBS. Standard:Ag12790. Range: 0.195-200 ng/mL



Cytometric bead array standard curve of MP50760-3, NTF3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60520-4-PBS. Detection antibody: 60520-1-PBS. Standard:Ag12790. Range: 1.563-200 ng/mL