

For Research Use Only

ND1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60534-1-PBS



Basic Information

Catalog Number: 60534-1-PBS	GenBank Accession Number: NC_012920	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 4535	CloneNo.: 2F8H6
Source: Mouse	UNIPROT ID: P03886	
Isotype: IgG1	Full Name: NADH dehydrogenase, subunit 1 (complex I)	
Immunogen Catalog Number: AG33365	Calculated MW: 36KD	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50754-1, ND1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60534-1-PBS. Detection antibody: 60534-2-PBS. Standard:Ag33365. Range: 0.391-100 ng/mL.