

For Research Use Only

# PDGFRB Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60535-5-PBS



## Basic Information

Catalog Number:

60535-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC032224

GeneID (NCBI):

5159

UNIPROT ID:

P09619

Full Name:

platelet-derived growth factor  
receptor, beta polypeptide

Calculated MW:

1106 aa, 124 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2F5F4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50756-3, PDGFRB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60535-5-PBS. Detection antibody: 60535-6-PBS. Standard:Eg0147. Range: 0.098-100 ng/mL.