

For Research Use Only

MCART1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 60544-1-PBS

Basic Information

Catalog Number:

60544-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG31357

GenBank Accession Number:

BC008500

GeneID (NCBI):

92014

UNIPROT ID:

Q9H1U9

Full Name:

mitochondrial carrier triple repeat 1

Purification Method:

Protein G Magarose purification

CloneNo.:

2E8C5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50768-1, MCART1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60544-1-PBS. Detection antibody: 60544-2-PBS. Standard: Ag31357. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50768-3, MCART1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60544-1-PBS. Detection antibody: 60544-4-PBS. Standard: Ag31357. Range: 0.098-100 ng/mL.