

For Research Use Only

NTRK1 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60548-3-PBS



Basic Information

Catalog Number:

60548-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27985

GenBank Accession Number:

NM_001007792

GeneID (NCBI):

4914

UNIPROT ID:

P04629

Full Name:

neurotrophic tyrosine kinase,
receptor, type 1

Calculated MW:

88 kDa

Purification Method:

Protein G purification

CloneNo.:

3A2G4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51082-1, NTRK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60548-2-PBS. Detection antibody: 60548-3-PBS. Standard:Ag27985. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP51082-2, NTRK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60548-3-PBS. Detection antibody: 60548-4-PBS. Standard:Ag27985. Range: 1.563-200 ng/mL