

For Research Use Only

CYP17A1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60554-1-PBS



Basic Information

Catalog Number:

60554-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5998

GenBank Accession Number:

BC062997

GeneID (NCBI):

1586

UNIPROT ID:

P05093

Full Name:

cytochrome P450, family 17,
subfamily A, polypeptide 1

Calculated MW:

57 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A5F2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50779-1, CYP17A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60554-1-PBS. Detection antibody: 60554-2-PBS. Standard:Ag5998. Range: 0.195-6.25 ng/mL